

# MobiLab®



“We previously had one or two labeling errors per month. Since going live with MobiLab®, there have been zero.”

**Deborah Murphy**

*Laboratory Information Systems Manager  
Hallmark Health System*

MobiLab® improves patient safety at your hospital. MobiLab barcode specimen collection brings positive patient identification, positive specimen identification, and specimen label printing to the bedside.



Smartphone



Zebra ZQ610  
Portable Printer



Mobile Carts and  
Bedside PCs



## Improve patient safety with mobile phlebotomy

By utilizing wireless handheld and printer technology, MobiLab is designed to reduce patient identification and specimen labeling errors, while increasing the productivity of your phlebotomy and nursing staff.



## Positive Patient ID

Barcode technology is used to positively identify a patient by wristband barcode at bedside. Phlebotomists also are required to confirm a secondary identification method in order to complete the Positive Patient ID process.



## Bedside Label Printing

Specimen collection labels for each patient are printed at the bedside through the use of portable label printers, eliminating preprinted label batches, a major source of mislabeled specimens.



## Positive Specimen ID

After tubes are labeled they are scanned to confirm the specimen label applied matches the patient scanned. Required collection time and collector ID are electronically documented. This simple step closes the loop on each collection event and ensures accurate documentation of collection details.

## Key Advantages Include:

- Eliminates patient and specimen identification errors
- Improves stat turnaround time
- Decreases unnecessary venipunctures via “Smart Draw” warnings that display previously collected specimens and future orders for the patient
- Provides the proper fit for each caregiver via multiple hardware options — handheld, laptop, PC, mobile cart — enabling a smooth transition and compliance
- Ranked #1 in the KLAS “Specimen Collection Barcoding” segment as published in the Top 20 Best in KLAS Awards report every year the segment existed: 2007, 2008, 2009, 2010, 2011, and 2012. ([www.KLASresearch.com](http://www.KLASresearch.com))



For more information about iatricSystems MobiLab or any other iatricSystems products or services, please contact using the information below.

\* Additional types of costs may include server and storage hardware, Microsoft licensing (OS, database, etc.), 3rd party licensing (digital certificates, backups, etc.), and 3rd party interface/ integration. Please consult with iatricSystems for the specific server hardware and software requirements for this product.